

Surgeon		Gilmer		Procedure		ACL Recon Quad Tendon Autograft Level 3	
Glove Size		8 Bio, 7 ½ white x 2		Position of Patient		Supine w/ lateral post. Patient's foot approx. 1 inch from end of bed. Unaffected leg wrapped in blanket and placed in middle of OR table	
Skin Prep		Chloraprep, Hibiclens & ETOH in Pacu, cover with blue towel. Clip hair in PACU pre-op					
Sutures Needles		ABS Tightrope (closed) AR1588TN BTB Tightrope AR1588BTB x1 ABS Tight Rope Button AR1588TB x1 Fiber Stick AR-7209 x1 Tiger Stick AR-7209T x1 #2 Fiberloop x2 3-0 Nylon FS-1 x2 0 Vicryl CT-1 pop-off 2-0 Vicryl CT-1		Drapes		Arthrotomy Pack, Extremity Sheet, Impervious U drapes x2, Impervious Stockinette x2, Coban, 6" Esmark, Mayostand cover, Suction Tubing, Mini C-Arm Drape, Double Basin set	
Dressings		¼" Steri Strips, Small Xeroform, 4x4s, 4" Webril, ABD, Double 6" Ace, Hinged ACL Brace		Instruments and Equipment			
Extras							
<p>Beginning: Marcaine 0.5% with Epi End: Ropivacaine 0.5% Plain 10ML + Morphine 5mg</p> <p>ANCEF (FOR ALLOGRAFT ONLY) Radiology: Have mini C-arm set up at beginning of case with pt. info. Takes two shots at end of case. Comes in same side</p> <p>Put shaver pedal on lift and Vulcan pedal on floor Bovie 40/40</p>				<p>Bone tray Stryker System 6 and one battery Arthroscopy Tray Arthrex Shaver Arthroscopy lens Crall ACL Set ACL work station Camera and light cord Pump and Patient tubing Medium Bump Marking pen 3.5 Arthrex Drill bit (knee cart, bottom shelf) Apollo RF 90 4.0 Dbl cut Cvd Quad Tendon Graft Cutting Blade, 10mm AR-2385-10 Sm. Black Triangle 10 blade</p>			
Notes				08/29/17 SB			
<p>Dr.Gilmer uses Mayo stand for graft station, not the mini table</p> <p>Quad AUTOGRAFT: Have #15 blade, Army Navy, Senns, Metz, Tenotomy, Allis, Arthrex quad harvester and Quad blade 10mm ready before scope. BTB tightrope, ABS Tightrope, ABS Button, Fiberloop x 2, Cut ruler at 2 cm and again at 9 cm. Work station: Graft sizer, 0 fiberwire x2, Needle Holder, Allis clamp, Adsons w/ teeth, Metz scissors, Marking pen, Ruler</p> <p>6 pictures per page, 1 copy and import to PACS</p>							